1.00579.0500
 1.07960.0500
 1.07961.0100
 1.07961.0500
 1.00869.0500
 1.01691.0025
 1.01691.0100
 1.03973.0001
 1.09016.0100
 1.09016.0500

Microscopy

DPX new

non-aqueous mounting medium for microscopy

Entellan™

rapid mounting medium for microscopy

Entellan™ new

rapid mounting medium for microscopy

Entellan™ new for cover slipper for microscopy

Canada balsam

for microscopy

M-GLAS®

liquid cover glass for microscopy

Neo-Mount™

anhydrous mounting medium for microscopy

For professional use only

IVD In Vitro Diagnostic Medical Device

Intended purpose

These ready-to-use anhydrous mounting media are suited for mounting dehydrated sample material of human origin after these have been fixed and embedded as necessary, and then histologically, bacteriologically, hematologically (enzyme-cytochemically), or cytologically stained and, where applicable, counterstained with other *in vitro* diagnostic products from our portfolio, thus making them evaluable for further diagnostic procedures. Samples are mounted on slides to enable the specimen material to be examined by light microscopy, at the same time preserving it and thus enabling it to be re-examined many years later.

The appropriate anhydrous mounting medium for the respective application is given in the corresponding instructions for use for our *In Vitro* Diagnostic staining solutions, solid dyes, and test kits.

Using the auxiliary reagents from our portfolio creates the conditions that enable authorized and qualified investigators to make a correct diagnosis at the end of the diagnostic process. In this regard, auxiliary IVD reagents serve *inter alia* to process human specimen material (e.g. fixing, decalcifying, dehydrating, clarifying, paraffin-embedding, mounting, microscoping, archiving). When used together with the corresponding staining solutions, this enables the visualization of cellular structures that are otherwise low in contrast, thus rendering them evaluable under the optical microscope. Further examinations may be necessary to reach a definitive diagnosis.

Principle

Mounting media are viscous, clear liquids with brilliant light-refraction properties. They are either won from natural materials or else are made e.g. of acryl-resin mixtures that are dissolved in aromatic solvents such as toluene, xylene, or a xylene substitute (e.g. Neo-Clear™, Cat. No. 1.09843).

In the last steps of the staining process prior to mounting, the still aqueous, stained specimen slides pass through a series of baths with ascending alcohol concentrations, ultimately ending up in an anhydrous solvent that is also referred to as an intermedium, e.g. toluene, xylene, or a xylene substitute (e.g. Neo-Clear™, Cat. No. 1.09843).

The anhydrous mounting media in their dissolved form are then dropped onto the stained and dehydrated specimen of human origin, and the slide is covered air-tight with a cover glass. The evaporation of the intermedium causes the mounting medium to harden, forming a solid, clear film under the cover glass, preserving the stained specimen material and thus enabling it to be kept for several years for re-analysis at a later date. As a result of the glass-similar refractive properties of the cover glass, the sample can now be observed under a microscope without any interference.

Thanks to the practical user-friendly dropping bottle, the mounting medium can be easily and safely dropped onto the slide without smearing. The closure of the nozzle ensures that the viscosity of the medium remains constant, meaning that the mounting medium is immediately ready for use.

Sample material

Starting materials are

- formalin-fixed, paraffin-embedded, histologically stained tissue specimens (3 - 5 µm thick paraffin sections)
- fixed and stained cytological smears, e.g. sputum, fine needle aspiration biopsies (FNAB), rinses, imprints, effusions
- air-dried, heat-fixed, and stained smears of bacteriological specimen material, e.g. liquid and solid enrichment cultures of bacteria from body fluids, exudates, pus
- hematologically processed and stained blood or bone-marrow smears from all regions of the human body.

Reagents

Cat. No. 1.00579	DPX new non-aqueous mounting medium for microscopy	500 ml
Cat. No. 1.07960	Entellan™ rapid mounting medium for microscopy	500 ml
Cat. No. 1.07961	Entellan™ new rapid mounting medium for microscopy	100 ml, 500 ml, 1 l
Cat. No. 1.00869	Entellan™ new for cover slipper for microscopy	500 ml
Cat. No. 1.01691	Canada balsam for microscopy	25 ml, 100 ml
Cat. No. 1.03973	M-GLAS® liquid cover glass for microscopy	500 ml
Cat. No. 1.09016	Neo-Mount™ anhydrous mounting medium for microscopy	100-ml dropping bottle, 500 ml

Specifications

Cat. No. 1.00579 - DPX new, non-aqueous mounting medium for microscopy

is a water-free mounting medium for microscopy, in which the teratogenic ingredient Dibutyl phthalate (DBP) has been avoided.

Refractive index (20°C)	1.518 - 1.521
Viscosity (20°C)	600 - 700 mPa*s

Cat. No. 1.07960 - Entellan™, rapid mounting medium for microscopy

is a water-free mounting medium for microscopy for the permanent mounting and storage of specimens, and consists of a polymer made of mixed acrylates dissolved in toluene. As it contains toluene, it should be used with water-free specimens that have been processed with xylene previous mounting.

Refractive index (20°C)	1.492 - 1.500
Density (20°C / 4°C)	0.925 - 0.935 g/cm ³
Viscosity (20°C)	60 - 100 mPa*s
Fluorescence	≤ 100 ppb

Cat. No. 1.07961 - Entellan™ new, rapid mounting medium for microscopy

is a water-free mounting medium for microscopy that consists of a polymer of mixed acrylates which are solubilized in xylene. Therefore, it should be used with specimens that have been cleared with xylene previous mounting.

Refractive index (20°C)	1.490 - 1.500
Density (20°C / 4°C)	0.94 - 0.96 g/cm ³
Viscosity (20°C)	250 - 600 mPa*s

Cat. No. 1.00869 - Entellan™ new for cover slipper for microscopy is a mounting medium for microscopy that is especially suited for standard commercial automated-mounting instruments that operate with glass coverslips. It is used as described in the instruction manual for cover slippers and the ideal amount of mounting agent is determined in a pilot run. There, empty cover glasses and specimen holders, according to the size of the cover glass and the size and thickness of the specimen, are used and these conditions are re-checked when a new bottle of the mounting medium is used. As its viscosity range is adjusted to a narrow range, the effort for new the calibration of the instrument is minimized.

Refractive index (20°C)	1.490 - 1.500
Viscosity (20°C)	500 - 600 mPa*s

Cat. No. 1.01691 - Canada balsam for microscopy

is a commonly used mounting medium for microscopy to prepare permanent slides. It is produced from the resin of the balsam fir tree and its use can be combined with xylene-containing specimens.

Refractive index (20°C)	1.515 - 1.530
Density (20°C / 4°C)	0.980 g/cm ³
Viscosity (20°C)	3000 mPa*s

Cat. No. 1.03973 - M-GLAS®, liquid cover glass for microscopy

is used in cytology instead of a cover glass to ensure that the stained specimens are homogeneously covered. A few drops are applied onto the specimen, taking care that the mounting medium is evenly distributed over the specimen material. After the solvent has evaporated, a solid, protective lacquer film remains that ensures that the specimen material is preserved. The M-GLAS® layer is not resilient to immersion oils. In exceptional cases, the time for which the specimen is exposed to the immersion oil should be kept to less than 10 minutes, since otherwise a residue-free removal of the oil can no longer be guaranteed. If the exposure time is longer, it is recommended to remove as much of the immersion oil from the specimen as possible, immerse it in xylene, and to mount it anew.

Refractive index (20°C)	1.490 - 1.500
Density (20°C / 4°C)	0.980 g/cm ³
Viscosity (20°C)	500 - 600 mPa*s
Fluorescence	≤ 250 ppb

Cat. No. 1.09016 - Neo-Mount™, anhydrous mounting medium for microscopy

is an extremely color-stable mounting medium for microscopy, which is produced with solvents based on mixtures of aliphatic hydrocarbons. It contains an aromatic-free substitute for xylene, thus, Neo-Mount™ needs to be combined with Neo-Clear™ (Cat. No. 1.09843) exclusively. Xylene must be avoided in the mounting step, as it will cause the slides to become cloudy and streaked. The application of Neo-Mount™ is not recommended in fluorescence microscopy for clinical diagnostics. In addition, by placing the dehydrated slides on filter paper for approx. 1 minute prior mounting, any excess of Neo-Clear™ could be circumvented, as air bubbles might arise under the coverslip. The same precondition should also be met when mounting specimens using cover-slip machines; in this area, Neo-Clear™ can be most efficiently eliminated by incubation of the slides for one minute in an empty slide rack.

Refractive index (20°C)	1.417 - 1.465
Viscosity (20°C)	250 - 350 mPa*s

Also required:

Cat. No. 1.00974	Ethanol denatured with about 1 % methyl ethyl ketone for analysis EMSURE®	1 l, 2.5 l
Cat. No. 1.08298	Xylene (isomeric mixture) for histology	4 l
Cat. No. 1.09843	Neo-Clear™ (xylene substitute) for microscopy	5 l, 25 l

Sample preparation

The sampling must be performed by qualified personnel.

All samples must be treated using state-of-the-art technology. All samples must be clearly labeled.

Suitable instruments must be used for taking samples and their preparation. Follow the manufacturer's instructions for application / use.

The specimen material is processed, stained (and counterstained where applicable), and mounted according to the instructions for use of our *In Vitro* Diagnostic staining solutions, solid dyes, and test kits.

Histological and cytological specimens must be completely dehydrated before mounting. In the last stage, either xylene or a xylene substitute should be used to prevent the occurrence of turbidity due to aqueous solutions.

Reagent preparation

All listed anhydrous mounting media are ready-to-use, dilution of the mounting media is not necessary.

When exchanging one anhydrous mounting medium in a cover slipper for another, e. g. when switching from Entellan™ to Entellan™ new, it is absolutely imperative to rinse the entire injection system of the cover slipper with the solvent xylene before using the new mounting medium. Only then can the new mounting medium be used.

If this is not done, oil-drop-shaped artefacts will form on the slide.

Procedure

The mounting medium must contain the same solvent, used for the water-clearing procedure to obtain the optimal optical properties and transparency of the slides.

All mounting procedures should be carried out in a fume hood.

The mounting medium is applied to the horizontal slide, using a glass rod or else directly dropping approx. 0.2 ml of one of the listed mounting media from the dropping bottle. As soon as a homogeneous distribution of the solution is guaranteed, gently add a clean cover glass, so that the space between the slide and cover glass is filled without air bubbles with mounting medium. Allow this setup to dry and harden for about 20 - 30 min in a horizontal position.

When pre-treated in the correct manner, the color of the specimens remains stable.

The use of immersion oil is recommended for the analysis of stained slides with a microscopic magnification >40x.

Note

In the case of mounted specimens, the cover slips can be detached again by immersing them in xylene. Specimens that have been mounted with M-GLAS® (Cat. No. 1.03973) can also be treated in this manner.

Cat. No. 1.00579 - DPX new, non-aqueous mounting medium for microscopy

Solvent	Xylene
Immersion time	approx. 65 hours

Cat. No. 1.07960 - Entellan™, rapid mounting medium for microscopy

Solvent	Xylene
Immersion time	approx. 24 hours

Cat. No. 1.07961 - Entellan™ new, rapid mounting medium for microscopy

Solvent	Xylene
Immersion time	approx. 72 hours

Cat. No. 1.00869 - Entellan™ new for cover slipper for microscopy

Solvent	Xylene
Immersion time	approx. 72 hours

Cat. No. 1.01691 - Canada balsam for microscopy

Solvent	Xylene
Immersion time	approx. 51 hours

Cat. No. 1.03973 - M-GLAS®, liquid cover glass for microscopy

Solvent	Xylene
Immersion time	approx. 17 hours

Cat. No. 1.09016 - Neo-Mount™, anhydrous mounting medium for microscopy

Solvent	Xylene
Immersion time	approx. 24 hours

Result

The use of these anhydrous, ready-to-use mounting media results in completely airtight specimen slides, the structure and stain pattern of which remains preserved over the long term, enabling them to be microscopically re-analyzed at a later date.

Trouble-shooting

Oil-drop-shaped artefacts on the slide

- When exchanging one anhydrous mounting medium in a cover slipper for another, e. g. when switching from Entellan™ to Entellan™ new, it is absolutely imperative to rinse the entire injection system of the cover slipper with the solvent xylene before using the new mounting medium. Only then can the new mounting medium be used.

Turbidity of the slides

- As a measure to ensure that the specimen slides retain optimal optical properties and their transparency, in all cases a mounting medium must be used that is based on the solvent / intermediate used for the clarification process. The Neo-Mount™ mounting medium is, for example, not compatible with xylene and hence should be used only in combination with the intermediate Neo-Clear™.

No color stability over longer storage times

- A minimum quality of the solvents must be observed. Technical-grade solvents may have a relatively high water content, which may result in incomplete dehydration and hence in the stained specimen becoming turbid or decolorized.
- Care must be taken to maintain a minimum quality and dye concentration of the staining solutions as a measure to stabilize the stain of the specimen.

Air bubbles and inclusions

- In all cases a mounting medium must be used that is based on the solvent / intermediate used for the clarification process.
- The volume of the mounting medium applied to the specimen must be carefully monitored (to avoid too much or too little mounting medium).
- The drying times for the specimens must be observed. The specimens must be completely dehydrated before microscopy with immersion oil, i. e. always allow the specimens to dry completely and mount thoroughly.
- The evaporation of the solvent after mounting must be borne in mind, and specimen slides must be dried for at least 20 - 30 min.

Technical notes

The microscope used should meet the requirements of a medical diagnostic laboratory.

When using automated-mounting instruments, please follow the instructions for use supplied by the supplier of the system and software.

Remove surplus immersion oil before filing.

Analytical performance characteristics

The present auxiliary reagents "DPX new", "Entellan™", "Entellan™ new", "Entellan™ new for cover slipper", "Canada balsam", "M-GLAS®", and "Neo-Mount™" aid in the microscopic examination of biological structures as described in the "Intended purpose" of this IFU. The use of the products is only to be carried out by authorized and qualified persons, this includes, among other things, sample and reagent preparation, sample handling, histoprocessing (Cat. Nos. 1.00579, 1.07960, 1.07961, 1.00869, 1.09016), decisions regarding suitable controls and more.

The analytical performance of the products is confirmed by testing each production batch.

For the following stains, the analytical performance was confirmed in terms of specificity, sensitivity and repeatability of the product with a rate of 100 %:

Cat. No. 1.00579 - DPX new

	Inter-assay Specificity	Inter-assay Sensitivity	Intra-assay Specificity	Intra-assay Sensitivity
Hematological, histological and physical methods				
Suitability for microscopy	10/10	10/10	7/7	7/7
Refractive index (n_{20}^D)	10/10	10/10	7/7	7/7

Analytical performance results

Cat. No. 1.07960 - Entellan™

	Inter-assay Specificity	Inter-assay Sensitivity	Intra-assay Specificity	Intra-assay Sensitivity
Hematological, histological and physical methods				
Suitability for microscopy	7/7	7/7	6/6	6/6
Refractive index (n_{20}^D)	7/7	7/7	6/6	6/6

Analytical performance results

Cat. No. 1.07961 - Entellan™ new

	Inter-assay Specificity	Inter-assay Sensitivity	Intra-assay Specificity	Intra-assay Sensitivity
Hematological, histological and physical methods				
Suitability for microscopy	20/20	20/20	6/6	6/6
Refractive index (n_{20}^D)	20/20	20/20	6/6	6/6

Analytical performance results

Cat. No. 1.00869 - Entellan™ new for cover slipper

	Inter-assay Specificity	Inter-assay Sensitivity	Intra-assay Specificity	Intra-assay Sensitivity
Hematological, histological and physical methods				
Suitability for microscopy	11/11	11/11	6/6	6/6
Refractive index (n_{20}^D)	11/11	11/11	6/6	6/6

Analytical performance results

Cat. No. 1.01691 - Canada balsam

	Inter-assay Specificity	Inter-assay Sensitivity	Intra-assay Specificity	Intra-assay Sensitivity
Hematological, histological and physical methods				
Suitability for microscopy	20/20	20/20	7/7	7/7
Refractive index (n_{20}^D)	20/20	20/20	7/7	7/7

Analytical performance results

Cat. No. 1.03973 - M-GLAS®

	Inter-assay Specificity	Inter-assay Sensitivity	Intra-assay Specificity	Intra-assay Sensitivity
Hematological, histological and physical methods				
Suitability for microscopy	9/9	9/9	7/7	7/7
Refractive index (n_{20}^D)	9/9	9/9	7/7	7/7

Analytical performance results

Cat. No. 1.09016 - Neo-Mount™

	Inter-assay Specificity	Inter-assay Sensitivity	Intra-assay Specificity	Intra-assay Sensitivity
Hematological, histological and physical methods				
Suitability for microscopy	9/9	9/9	7/7	7/7
Refractive index (n_{20}^D)	9/9	9/9	7/7	7/7

Analytical performance results

Intra- (performed on the same batch) and inter-assay (performed on different batches) data list the number of correctly stained structures in relation to the number of performed assays.

The results of this Performance Evaluation confirms that the products are suitable for the intended use and performs reliably.

Diagnostics

Diagnoses are to be made only by authorized and qualified personnel. Valid nomenclatures must be used.

These products are auxiliary reagents these, when used together with other IVD products such as staining solutions, render human specimen material evaluable for diagnostic purposes.

Further tests must be selected and implemented according to recognized methods.

Suitable controls should be conducted with each application in order to avoid an incorrect result.

Storage

Cat. Nos. 1.00579, 1.07960, 1.07961, 1.00869, 1.03973, 1.09016: Store the listed mounting media at +15 °C to +25 °C.

Cat. No. 1.01691:

Store the Canada balsam for microscopy at +5 °C to +30 °C.

Shelf-life

Cat. Nos. 1.00579, 1.07960, 1.07961, 1.00869, 1.03973, 1.09016: The listed mounting media can be used until the stated expiry date.

After first opening of the bottle, the contents can be used up to the stated expiry date when stored at +15 °C to +25 °C.

Cat. No. 1.01691:

The Canada balsam for microscopy can be used until the stated expiry date. After opening the bottle the first time and subsequent storage of the tightly reclosed bottle at +5°C to +30 °C, the medium can be used up to the printed expiry date.

Additional instructions

For professional use only.

In order to avoid errors, the application must be carried out by qualified personnel only.

National guidelines for work safety and quality assurance must be followed. Microscopes equipped according to the standard must be used.

Protection against infection

Effective measures must be taken to protect against infection in line with laboratory guidelines.

Instructions for disposal

The package must be disposed of in accordance with the current disposal guidelines.

Used solutions and solutions that are past their shelf-life must be disposed of as special waste in accordance with local guidelines. Information on disposal can be obtained under the Quick Link "Hints for Disposal of Microscopy Products" at www.microscopy-products.com. Within the EU the currently applicable REGULATION (EC) No 1272/2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006 applies.

Auxiliary reagents

Cat. No. 1.00974 Ethanol denatured with about 1 l, 2.5 l
1 % methyl ethyl ketone
for analysis EMSURE®

Cat. No. 1.00983 Ethanol absolute for analysis EMSURE® 1 l, 2.5 l, 5 l
ACS, ISO, Reag. Ph Eur

Cat. No. 1.04699	Immersion oil for microscopy	100-ml dropping bottle, 100 ml, 500 ml
Cat. No. 1.08298	Xylene (isomeric mixture) for histology	4 l
Cat. No. 1.09843	Neo-Clear™ (xylene substitute) for microscopy	5 l, 25 l

Hazard classification

Cat. Nos. 1.00579, 1.07960, 1.07961, 1.00869, 1.01691, 1.03973, 1.09016

Please observe the hazard classification printed on the label and the information given in the safety data sheet.

The safety data sheet is available on the website and on request. CAUTION! Cat. Nos. 1.00579, 1.07960, and 1.03973 contain CMR substances. Please observe the corresponding safety instructions given in the safety data sheet.

Main components of the products

Cat. No. 1.00579

Copolymer in 70 % (w/w) xylene

Cat. No. 1.07960

Mixed acrylate in 75 % (w/w) toluene

1 l = 0.93 kg

Cat. No. 1.07961

Polymer of mixed acrylates in 60 % (w/w) xylene

1 l = 0.95 kg

Cat. No. 1.00869

Polymer of mixed acrylates in 60 % (w/w) xylene

1 l = 0.95 kg

Cat. No. 1.01691

CAS-No 8007-47-4

1 l = 0.98 kg

Cat. No. 1.03973

Polymer of mixed acrylates in 73.3 % (w/w) toluene

1 l = 0.91 kg

Cat. No. 1.09016

Polymer of mixed acrylates in 64 % (w/w) ShellSol 140/165

Other IVD products

Cat. No. 1.00496	Formaldehyde solution 4%, buffered, pH 6.9 (approx. 10% Formalin solution) for histology	350 ml and 700 ml (in bottle with wide neck), 5 l, 10 l, 10 l Titripac®
Cat. No. 1.01646	PAS staining kit for detection of aldehyde and mucosubstances	2x 500 ml
Cat. No. 1.05174	Hematoxylin solution modified acc. to Gill III for microscopy	500 ml, 1 l, 2.5 l
Cat. No. 1.09204	Giemsa's azur eosin methylene blue solution for microscopy	100 ml, 500 ml, 1 l, 2.5 l
Cat. No. 1.11609	Histosec™ pastilles solidification point 56-58°C embedding agent for histology	1 kg, 10 kg (4x 2.5 kg), 25 kg
Cat. No. 1.11885	Gram-Color stain set for the Gram staining method	1 set
Cat. No. 1.15161	Histosec™ pastilles (without DMSO) solidification point 56-58°C embedding agent for histology	10 kg (4x 2.5 kg), 25 kg

General remark

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and / or its authorised representative and to your national authority.

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Cat. No. 1.00579

H225: Highly flammable liquid and vapor.

H304: May be fatal if swallowed and enters air-ways.

H312 + H332: Harmful in contact with skin or if inhaled.

H315: Causes skin irritation.

H319: Causes serious eye irritation.

H335: May cause respiratory irritation.

H412: Harmful to aquatic life with long lasting effects.

P210: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

P273: Avoid release to the environment.

P280: Wear protective gloves/ protective clothing/ eye protection/ face protection.

P301 + P310: IF SWALLOWED: Immediately call a POISON CENTER/ doctor.

P303 + P361 + P353: IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water.

P331: Do NOT induce vomiting.

Cat. No. 1.07960

H225: Highly flammable liquid and vapor.

H304: May be fatal if swallowed and enters air-ways.

H315: Causes skin irritation.

H336: May cause drowsiness or dizziness

H361d: Suspected of damaging the unborn child.

H373: May cause damage to organs (Central nervous system) through prolonged or repeated exposure.

H412: Harmful to aquatic life with long lasting effects.

P202: Do not handle until all safety precautions have been read and understood.

P210: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

P273: Avoid release to the environment.

P301 + P310: IF SWALLOWED: Immediately call a POISON CENTER/ doctor.

P303 + P361 + P353: IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water.

P331: Do NOT induce vomiting.

Cat. No. 1.07961

H226: Flammable liquid and vapor.

H304: May be fatal if swallowed and enters air-ways.

H312 + H332: Harmful in contact with skin or if inhaled.

H315: Causes skin irritation.

H317: May cause an allergic skin reaction.

H319: Causes serious eye irritation.

H335: May cause respiratory irritation.

H373: May cause damage to organs (Central nervous system, Liver, Kidney) through prolonged or repeated exposure.

H412: Harmful to aquatic life with long lasting effects.

P210: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

P273: Avoid release to the environment.

P280: Wear protective gloves/ protective clothing/ eye protection/ face protection/ hearing protection.

P301 + P310: IF SWALLOWED: Immediately call a POISON CENTER/ doctor.

P303 + P361 + P353: IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water.

P331: Do NOT induce vomiting.

Cat. No. 1.00869

H226: Flammable liquid and vapor.
 H304: May be fatal if swallowed and enters air-ways.
 H312 + H332: Harmful in contact with skin or if inhaled.
 H315: Causes skin irritation.
 H319: Causes serious eye irritation.
 H335: May cause respiratory irritation.
 H373: May cause damage to organs (Central nervous system, Liver, Kidney) through prolonged or repeated exposure.
 H412: Harmful to aquatic life with long lasting effects.

P210: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

P273: Avoid release to the environment.

P280: Wear protective gloves/ protective clothing/ eye protection/ face protection/ hearing protection.

P301 + P310: IF SWALLOWED: Immediately call a POISON CENTER/ doctor.

P303 + P361 + P353: IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water.

P331: Do NOT induce vomiting.

Cat. No. 1.01691

H226: Flammable liquid and vapor.

P210: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

P233: Keep container tightly closed.

P240: Ground and bond container and receiving equipment.

P241: Use explosion-proof electrical/ ventilating/ lighting/ equipment.

P242: Use non-sparking tools.

P243: Take action to prevent static discharges.

Cat. No. 1.03973

H225: Highly flammable liquid and vapor.

H304: May be fatal if swallowed and enters airways.

H315: Causes skin irritation.

H336: May cause drowsiness or dizziness.

H361: Suspected of damaging fertility or the unborn child.

H373: May cause damage to organs (Central nervous system) through prolonged or repeated exposure.

H412: Harmful to aquatic life with long lasting effects.

P210: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

P233: Keep container tightly closed.

P273: Avoid release to the environment.

P301 + P310: IF SWALLOWED: Immediately call a POISON CENTER/ doctor.

P314: Get medical advice/ attention if you feel unwell.

P331: Do NOT induce vomiting.

Cat. No. 1.09016

H225: Highly flammable liquid and vapor.

H336: May cause drowsiness or dizziness.

H412: Harmful to aquatic life with long lasting effects.

P210: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

P233: Keep container tightly closed.

P240: Ground and bond container and receiving equipment.

P241: Use explosion-proof electrical/ ventilating/ lighting/ equipment.

P242: Use non-sparking tools.

P273: Avoid release to the environment.

EUH066: Repeated exposure may cause skin dryness or cracking.

Consult instructions
for use

Manufacturer

Catalog number

Batch code

Caution, consult
accompanying documents

Use by
YYYY-MM-DD

Temperature
limitation

Status: 2023-May-17

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

© 2023 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Merck and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly available resources.

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